

Technical Service Bulletin			
Edition: Final V1	TSB#: (ECR1871&1903)	Issue date: April 06, 2023	Page 1 of 2

Title: Peristaltic Pump & Support Standardization

Old parts no longer used / New parts to be implemented:

OLD Part – no longer available:	NEW Part – now available:	For Machine Types:
60000046 Peristaltic Pump 3 Rolls 24V	60000044 Peristaltic Pump 24V 4 Rolls	Diamant Pro Onyx Pro
60000049 Peristaltic Pump 24V		
260673_Peristaltic_Pump_Support	21004128_Chassis_Support_Peristaltic_Pump	

Applies to machine types:

- Diamant Pro
- Onyx Pro

Implemented from first serial number:

IMPLEMENTATION	Diamant Pro	Onyx Pro
New Peristaltic 60000044	To be implemented	To be implemented
New Support 21004128	To be implemented	To be implemented

Change information:

Change peristaltic & support to standardize auto steam solution between the different platforms.

Both changes are backwards compatible.

We recommend adjusting the pump speed in the technical menu after installing the new peristaltic pump, the settings are similar, but it is recommended to make a double check.

EX3 is already fitted with new parts since SmartSteam project launch.

Backward compatibility:

The new part(s) are fully backward compatible with the previous version machines.

Machines in the field:

For machines in the field, require no change.

Warranty information:

No warranty changes.

Spare part information:

The new part(s) will be the spare part(s) offered for all after sales services.

Technical Service Bulletin



Edition: Final V1

TSB#: (ECR1871&1903)

Issue date: April 06, 2023

Page 2 of 2

Improvement descriptions:

Old Version: 6000049	New Version: 6000044
	
<p data-bbox="347 1070 625 1102">Old Version: 6000049</p> 